


A Magnetizing Technology:

Nanowire Biosensors

Kamila Kunes and Ben Yang
Electrical Engineering
Group 2


Nanowire Biosensors So Far...



Labeled:

- Pre-Amplifying signal
- Indirect measurement


Label Free





- **Noise** from solution
- **Debye** length / charge screening
- **Fabrication** random

Aptamers require a special touch to develop

Aptamer Approach



But Wait... Why Should We Care?



An Alternative? → Magnetic Nanowire Biosensors


Fundamentally:

1. Target exists in media
2. Target changes resistance
3. Measure resistance change


How does target change resistance?

Nanowire FET → Electrically

Permalloy nanowire → Magnetically



Nanowire FET




Permalloy Nanowire

Why Magnetic Systems?

At its core: **INTERFERENCE**


- low background **noise**
- competitive **sensitivity**
- large dynamic **range**
- low **cost**
- capability of **multiplexed** detection




Magnetic Nanowire Biosensors – an overview

Vary σ with magnetic field effects:


- Giant Magnetoresistance (electron spin)
- Domain Wall Shift (notched nanowire)
- Magnetophoresis (zigzag nanowire)







Domain Wall Shift



Zigzag NW



MBS – Nanobeads as Labels



Giant magnetoresistive-based biosensing probe station system for multiplex protein assays
Yi Wang ¹ , Wei Wang ¹ , Lina Yu, Liang Tu, Yinglong Feng, Todd Klein, Jian-Ping Wang


Magnetic Particle Labeling

Permalloy (Ni₈₀Fe₂₀) nanowire

Search currents measure resistance changes in nanowires

0.3 T magnetic field +y
ensure that all the biosensor configurations are fully saturated

Allows for single label detection




$r = 500\text{nm}$ $L = 5\mu\text{m}$
 $W = 500\text{nm}$ $t = 20\text{nm}$

Labels are about 200nm wide


$$R(\theta) = R_{\text{perp}} + \delta R \cos^2(\theta)$$

(a)



R_{perp} = resistance of notch when the magnetization is perpendicular to the current ($\theta = 90^\circ$)


(b)



δR = change in resistance between 0° and 90°

Domain wall structure in larger cross sections influence resistance to a lesser extent




(c)




Closed domain structure reduces probability of label-label interactions




Zigzag Magnetic Label for Detection of Cancer Cells

Sensor size	Magnetic label particle(s)
150–800 nm nanowire	Fe_3O_4
Label size (diameter)	Normalized signal magnitude ^a
10 nm	1.74×10^{13}





Magnetic Nanowire Labels



Future - Magnetic Systems for Streamlined Surface



1. Diagnostics
 - a. Magnetic labeling for magnetic nanowire sensors
2. Treatment - Hyperthermia
 - a. Vibrate nanowires within cancer cells to “cook” them
3. Checkups

Thank You!

Article Links

- <https://www.intechopen.com/chapters/57318> 2017
- <https://pubs.rsc.org/en/content/articlelanding/2013/lc/c3lc50457c> 2013
- <https://iopscience.iop.org/article/10.1088/1361-6528/aba1ba> 2020
- <https://conservancy.umn.edu/handle/11299/226637> 2021
- <https://experts.umn.edu/en/publications/selective-detection-of-cancer-cells-using-magnetic-nanowires> 2021
- <https://ieeexplore.ieee.org/document/6655995> 2014